


Distal Volar Radial Locking Plate System

Product Information

Plates


Distal Volar Radial Plate, with drill guide

REF. NO.	Direction	Holes
026741003	Large, Left	3
026742003	Large, Right	3
026741004	Large, Left	4
026742004	Large, Right	4
026741005	Large, Left	5
026742005	Large, Right	5
026741007	Large, Left	7
026742007	Large, Right	7
026741009	Large, Left	9
026742009	Large, Right	9
026751003	Small, Left	3
026752003	Small, Right	3
026751004	Small, Left	4
026752004	Small, Right	4
026751005	Small, Left	5
026752005	Small, Right	5
026751007	Small, Left	7
026752007	Small, Right	7
026751009	Small, Left	9
026752009	Small, Right	9



Distal Volar Radial Locking Plate, with drill guide, combi-holes

REF. NO.	Direction	Holes
026801003	Large, Left	3
026802003	Large, Right	3
026801004	Large, Left	4
026802004	Large, Right	4
026801005	Large, Left	5
026802005	Large, Right	5
026801007	Large, Left	7
026802007	Large, Right	7
026801009	Large, Left	9
026802009	Large, Right	9
026821003	Small, Left	3
026822003	Small, Right	3
026821004	Small, Left	4
026822004	Small, Right	4
026821005	Small, Left	5
026822005	Small, Right	5
026821007	Small, Left	7
026822007	Small, Right	7
026821009	Small, Left	9
026822009	Small, Right	9



Distal Volar Radial Locking Plate

REF. NO.	Direction	Holes
026071003	Left	3
026072003	Right	3
026071004	Left	4
026072004	Right	4
026071005	Left	5
026072005	Right	5
026071007	Left	7
026072007	Right	7
026071009	Left	9
026072009	Right	9



Distal Volar Radial Locking Plate System



DOUBLE MEDICAL TECHNOLOGY INC.

No. 18, Shanjianhong East Road, Haicang District,
 Xiamen, 361022, P.R. China.
 Tel: +86 592 6087101 Fax: +86 592 6587078
 E-mail: info@double-medical.com
 Website: www.double-medical.com



Watershed Line

Contoured to match the watershed line, the distal end of the plate can provide a visual guide for optimal placement and is low-profile to blend into the bone which may mitigate the risk of tendon irritation.



Easy distinguish between L/R:

Left with golden sleeve
Right with blue sleeve

FAST Guide Sleeve

Technology allows for easy drilling of fixed angle locking screws and visually distinguishes left and right plates.

Narrower Shaft

Increases the ease of fitting the plate to the bone while still providing more fixation options than the current other anatomical plate.



Low profile construct

Bevelled distal edge, rounded and polished surface and countersunk screws for reduced soft tissue irritation.

Bendable tabs

If necessary, tabs can be bent in order to suit the individual anatomical conditions of the bone.

Additional distal screw options

Enable support of radial styloid, lunate facet and DRUJ.



Warning

This flyer is just for understanding the specific product features. For clinical usage, please refer to the surgical guide. Instruction by experienced surgeon is highly recommended.

Indications

- The complex intra-articular and extra-articular fractures of the distal radius.